

TERMINAL URINE COLLECTION **(for Schistosomiasis detection)**

Obtain a sterile 50 ml yellow top container from your doctor or pathology collection centre

Please note

- A urine sample for Schistosomiasis detection can be collected from any urine you pass at any time of the day.

Instructions

1. Write your full name and date of birth on the container
2. Pass most of the urine into the toilet, and then stop the flow of urine.
3. Collect the final part of the urine you pass into the container provided.
4. Secure the lid.
5. Write the date and time of the collection on the container
6. Return the container to your doctor or to the pathology collection centre (with a referral form).

Note:

The specimen can be refrigerated if there is any delay